

Finnish PHR Calories Burned profile:

Name		Card.	Type	Remarks
extension				
	creatingApplication	0..1		
code			CodeableConcept	Binding: Value Set Finnish PHR Fitness Results (required)
	coding			
	caloriesBurned	1..1		
	--version	0..0		disabloitu
	--userSelected	0..0		disabloitu
	caloriesBurned24h	0..1		
	--version	0..0		disabloitu
	--userSelected	0..0		disabloitu
	weeklycaloriesBurned	0..1		
	--version	0..0		disabloitu
	--userSelected	0..0		disabloitu
valueQuantity				
	comparator			disabloitu
bodySite				disabloitu
method				disabloitu

Finnish PHR Distance profile:

Name		Card.	Type	Remarks
extension				
	creatingApplication	0..1		
code			CodeableConcept	Binding: Value Set Finnish PHR Fitness Results (required)
	coding			
	distance	1..1		
	--version	0..0		disabloitu
	--userSelected	0..0		disabloitu
	distance24h	0..1		
	--version	0..0		disabloitu
	--userSelected	0..0		disabloitu
	weeklyDistance	0..1		
	--version	0..0		disabloitu
	--userSelected	0..0		disabloitu
valueQuantity				
	comparator			disabloitu
bodySite				disabloitu
method				disabloitu

Finnish PHR Moderate to Vigorous Physical Activity Time profile:

Name		Card.	Type	Remarks
extension				
	creatingApplication	0..1		
code			CodeableConcept	Binding: Value Set Finnish PHR Fitness Results (required)
	coding			
	moderateToVigorousPAweek	0..1		
	--version	0..0		disabloitu
	--userSelected	0..0		disabloitu
	moderateToVigorousPA24h	0..1		
	--version	0..0		disabloitu
	--userSelected	0..0		disabloitu
valueQuantity				
	comparator	0..1		kuvattava paremmin
bodySite				disabloitu
method				disabloitu

Finnish PHR Steps taken profile:

Name		Card.	Type	Remarks
extension				
	creatingApplication	0..1		
code			CodeableConcept	Binding: Value Set Finnish PHR Fitness Results (required)
	coding			
	steps	1..1		
	--version	0..0		disabloitu
	--userSelected	0..0		disabloitu
	dailySteps	0..1		
	--version	0..0		disabloitu
	--userSelected	0..0		disabloitu
	weeklySteps	0..1		
	--version	0..0		disabloitu
	--userSelected	0..0		disabloitu
valueQuantity				
	comparator			disabloitu
	unit			OK
bodySite				disabloitu
method				disabloitu